

Bringing bright ideas to light.



Exam, Procedure & Surgical Lighting

MI-750

The ideal mid-sized LED light for minor surgery, procedures and exams. The MI-750 delivers maximum brightness, sharpness, superb color accuracy and great articulation.


The best choice for brighter
Procedure outcomes.

3-stage dimming via the
sterilizable handle lets you
easily control brightness &
on/off function from within
the sterile field.



MI-750

Bringing bright ideas to light.



MI-750 Specifications

Optical

| Parameter | MI-750 |
|-----------------------|-------------------------|
| # of LED's | 18 |
| Output | 75,000 Lux @ 39.4" (1m) |
| CCT | CCT 4,300°K |
| CRI | 97 |
| Focal Length | 39.4" (1m) |
| Dimming Levels | 3 |
| Dimming Percentages | 100%, 75%, 50% |
| Endo Mode | No |
| Light Field (d10) | ≈Ø7.7" (Ø196mm) |
| Depth of Illumination | 26.5" (673mm) |
| Irradiance | 196 W/m ² |

Electrical Specifications

| | |
|----------------|-----------------|
| Supply Circuit | 100-240VAC, 65W |
|----------------|-----------------|

Mechanical Specifications

| | |
|--------------------------------------|--|
| Light Head Dimensions | Ø14.5" (368mm) 3.25" H (83.3mm) |
| Arm Assembly Weight | Single: 48.3 lbs (22.0 kg) Dual: 77.1 lbs (35.0 kg) |
| Articulating Arm Vertical Movement | +15° Up, -60° Down |
| Articulating Arm Horizontal Movement | Movement ≈540° |

A mounting configuration for every setting



Floor



Dual Ceiling Mount



Ceiling Mount

Getting a better handle on things.

Controlling the brightness and position of your lights while maintaining sterility is critical. With the patented design of the MI Series lights - it's incredibly easy. Just tap the button on the bottom of the sterilizable handle to adjust brightness or turn the light on & off.



3 year warranty

Uncompromising quality, superb ergonomics.

| Order Code | Description |
|------------|-----------------------------------|
| 061422 | MI-750 Wall-Mount style |
| 061423 | MI-750 Mobile Floor-Mount style |
| 061424* | MI-750 Single Ceiling Mount style |
| 061425* | MI-750 Dual Ceiling Mount style |

* Standard mounting hardware kits accommodate ceilings up to 9'-8" high. Downrods for ceiling heights up to 12'-2" are available at extra charge

Call or click for ordering and information

800-862-7616
valleymed.ca

