

Bringing bright ideas to light.



LED & Fiberoptic Headlights & Lightsources

CHROME MC7 Pro

Surgical LED Headlight System

Comfortable, Light and Bright. Completely cordless and amazingly bright, the MC7 Pro provides the illumination and long-life battery capability required for all-day surgical use.

Air Strap
Headband (ST)



Hard Top
Headband (HT)



Uncompromising quality,
superb ergonomics.



Shown in optional ratchet top style headband, (HT)

A choice of 2 headband styles is available with an easily positioned headlight module, and optional clip-on magnifiers / loupes for excellent visualization give surgeons the freedom, comfort and capability to concentrate on their cases.

Bringing bright ideas to light.

MC7 Pro

Comfortable/cordless/
bright white light

Battery Charge Indicator
Visible charge level LED indicators

Long Life Batteries
6.5 - 7 hours (Infinite with hot-swap)

Comfortable
New Slimstrap™ design, for ideal balance

Digital Dimming
Continuously adjustable from 0-100%

Interchangeable Post Design
allows multiple light modules

High Intensity
S/B 275,000 Lux adjustable spot Size 9 - 20 cm (3.6" - 7.8")

Optional Magnifiers

Quiet-Cool™

Quiet-Cool™ heat dispensing technology both reduces temperature of the LED and reduces bezel size, keeping the Chrome Pro™ both cool and refined.

Battery


Dual Lithium Polymer battery packs provide constant steady brightness and longevity to the headlight operation while maintaining Hot-Swap™ capabilities.

Charging

The Chrome Pro™ onboard unit charger gives you the versatility of charging the headlight from any USB power supply. Choose the optional 4-bay desktop charging station and some extra batteries to be assured of 'always ready' performance.

Headlight Module

A two-pivot arm design allows the headlight module to be easily positioned for optimal working position and beam angle.



Order Code	Description
VMCH-07-HT	CHROME MC7 PROTM O.R. LED Headlight W/Slimstrap Ratchet Top Adj.
VMCH-07-HT-HK	Hospital Kit - CHROME MC7 PROTM O.R. LED Headlight W/Slimstrap Ratchet Top Adj. Kit adds 4 extra batteries (Total 6) , 4-Bay desktop charger & 2D clip-on magnifiers
VMCH-07-ST	CHROME MC7 PROTM O.R. LED Headlight W/Air Strap (Velcro) Top Adj.
VMCH-07w-ST-HK	Hospital Kit - CHROME MC7 PROTM O.R. LED Headlight W/Air Strap (Velcro) Top Adj. Kit adds 4 extra batteries (Total 6) , 4-Bay desktop charger & 2D clip-on magnifiers
Order Code	Accessories
VM-LIPO-2C	Batteries - Set of 2 - Chrome series
VMS01-F	Amber UV filter
VMDT-04	4-bay desktop charger
VMM01-MXX	Magnifying loupes - (2, 2.5, 3, 3.5 4, 4.5diopter)

Replacement padding sets available in Neoprene Velcro & Soft Pillow Velcro

Call or click for ordering and information

800-862-7616
valleymed.ca

